

TransNGS[®] Library Amplification SuperMix

Cat. No. KA101

Storage: at -20°C for two years

Description

TransNGS[®] Library Amplification SuperMix is optimized for high fidelity, robust amplification of next generation sequencing (NGS) libraries regardless of GC contents. The ready-to-use 2×Mix contains DNA polymerase, dNTPs and optimized buffer, and requires only the addition of primers and template.

Highlights

- High fidelity amplification.
- Low amplification bias.
- High sensitivity and high specificity.
- Hot start

Applications

- Next-generation sequencing library amplification.

Kit Contents

Component	KA101-01	KA101-02
TransNGS [®] Library Amplification SuperMix	1 ml	5×1 ml
Nuclease-free Water	1 ml	5 ml

Reaction Components

Component	Volume	Final Concentration
Adapter-ligated DNA	0.5-4 µl	-
Library Amplification Forward Primer (10 µM)	0.4 µl	0.2 µM
Library Amplification Reverse Primer (10 µM)	0.4 µl	0.2 µM
TransNGS [®] Library Amplification SuperMix	10 µl	1×
Nuclease-free Water	Variable	-
Total volume	20 µl	-

Recommended thermal cycling conditions

98°C	3 min	} 2-15 cycles**
98°C	30 sec	
x°C*	30 sec	
72°C	30 sec	
72°C	3 min	
≤10°C	Hold	

*Depending on the PCR primer length and GC contents.

**Depending on the amount of the starting material and your library preparation method.



Notes

- All components should be thawed and mixed thoroughly before use.
- We suggest to purify DNA after adapter ligation. Higher yield will be obtained with high quality DNA template.
- DNA polymerase cannot incorporate dUTP. dUTP-containing primers or templates in the reaction are not recommended.

